

## Anti-TGM1 antibody

Cat. No.	ml161172
Package	25 μl/100 μl/200 μl
Storage	-20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview		
Description	Anti-TGM1 rabbit polyclonal antibody	
Applications	ELISA, IHC	
Immunogen	Synthetic peptide of human TGM1	
Reactivity	Human	
Content	0.2 mg/ml	
Host species	Rabbit	
lg class	Immunogen-specific rabbit IgG	
Purification	Antigen affinity purification	
Target information		
Symbol	TGM1	
Full name	transglutaminase 1 (K polypeptide epidermal type I, protein-glutamine-gamma- glutamyltransferase)	
Synonyms	LI; KTG; LI1; TGK; ICR2; ARCI1; TGASE	
Swissprot	P22735	
Target Background		

## **Target Background**


The protein encoded by this gene is a membrane protein that catalyzes the addition of an alkyl group from an akylamine to a glutamine residue of a protein, forming an alkylglutamine in the protein. This protein alkylation leads to crosslinking of proteins and catenation of polyamines to proteins. This gene contains either one or two copies of a 22 nt repeat unit in its 3' UTR. Mutations in this gene have been associated with autosomal recessive lamellar ichthyosis (LI) and nonbullous congenital ichthyosiform erythroderma (NCIE).



订购热线: 4008-898-798

Applications Immunohistochemistry

Predicted cell location: Membrane Positive control: Human cervical cancer Recommended dilution: 25-100



The image on the left is immunohistochemistry of paraffin-embedded Human cervical cancer tissue using ml161172(TGM1 Antibody) at dilution 1/25, on the right is treated with synthetic peptide. (Original magnification: ×200)

## ELISA

Recommended dilution: 2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695,2881505696

邮箱: mlbio\_cn@yeah.net 网址: www.mlbio.cn