

订购热线: 4008-898-798

Anti-APBA2 antibody

Cat. No.	ml161289
Package	25 μl/100 μl/200 μl
Storage	-20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Q99767

Product overview	
Description	Anti-APBA2 rabbit polyclonal antibody
Applications	ELISA, IHC
Immunogen	Synthetic peptide of human APBA2
Reactivity	Human
Content	0.6 mg/ml
Host species	Rabbit
lg class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification
Target information	
Symbol	APBA2
Full name	amyloid beta (A4) precursor protein-binding, family A, member 2
Synonyms	X11L; MINT2; LIN-10; HsT16821; X11-BETA; D15S1518E; MGC:14091

Target Background

Swissprot

The protein encoded by this gene is a member of the X11 protein family. It is a neuronal adapter protein that interacts with the Alzheimer's disease amyloid precursor protein (APP). It stabilizes APP and inhibits production of proteolytic APP fragments including the A beta peptide that is deposited in the brains of Alzheimer's disease patients. This gene product is believed to be involved in signal transduction processes. It is also regarded as a putative vesicular trafficking protein in the brain that can form a complex with the potential to couple synaptic vesicle exocytosis to neuronal cell adhesion. Multiple transcript variants encoding different isoforms have been found for this gene.


reke



订购热线: 4008-898-798

Applications Immunohistochemistry

Predicted cell location: Cytoplasm Positive control: Human liver cancer Recommended dilution: 100-300





The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using ml161289(APBA2 Antibody) at dilution 1/80, on the right is treated with synthetic peptide. (Original magnification: ×200)

ELISA

Recommended dilution: 2000-5000

Predicted cell location: Cytoplasm Positive control: Human thyroid cancer Recommended dilution: 100-300





The image on the left is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using ml161289(APBA2 Antibody) at dilution 1/80, on the right is treated with synthetic peptide. (Original magnification: ×200)

- 联系电话: 4008-898-798, 021-61725725
- 联系QQ: 2881505695,2881505696、
- 邮箱: mlbio_cn@yeah.net 网址: www.mlbio.cn