

## Anti-NOX3 antibody

 Cat. No.
 ml161971

 Package
 25 μl/100 μl/200 μl

 Storage
 -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview Description Applications Immunogen Reactivity Content Host species Ig class Purification

Anti-NOX3 rabbit polyclonal antibody ELISA, IHC Synthetic peptide of human NOX3 Human 0.8 mg/ml Rabbit Immunogen-specific rabbit IgG Antigen affinity purification

Target information Symbol Full name Synonyms Swissprot NOX3 NADPH oxidase 3 MOX-2; GP91-3 Q9HBY0

## **Target Background**

This gene encodes a member of the NOX family of NADPH oxidases. These enzymes have the capacity to generate superoxide and other reactive oxygen species (ROS) and transport electrons across the plasma membrane. The ROS generated by family members have been implicated in numerous biological functions including host defense, posttranlational processing of proteins, cellular signaling, regulation of gene expression, and cell differentiation.



订购热线: 4008-898-798

## Applications Immunohistochemistry

Predicted cell location: Cytoplasm Positive control: Human brain Recommended dilution: 50-200





The image on the left is immunohistochemistry of paraffin-embedded Human brain tissue using ml161971(NOX3 Antibody) at dilution 1/40, on the right is treated with synthetic peptide. (Original magnification:

×200)

## ELISA

Recommended dilution: 2000-5000

Predicted cell location: Cytoplasm Positive control: Human esophagus cancer Recommended dilution: 50-200





The image on the left is immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using ml161971(NOX3 Antibody) at dilution 1/40, on the right is treated with synthetic peptide. (Original magnification: ×200)

- 联系电话: 4008-898-798, 021-61725725
- 联系QQ: 2881505695,2881505696、
- 邮箱: mlbio\_cn@yeah.net 网址: www.mlbio.cn