

## Anti-CCDC85C antibody

|                 |                                                         |
|-----------------|---------------------------------------------------------|
| <b>Cat. No.</b> | ml163389                                                |
| <b>Package</b>  | 25 µl/100 µl/200 µl                                     |
| <b>Storage</b>  | -20°C, pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol |

### Product overview

|                     |                                         |
|---------------------|-----------------------------------------|
| <b>Description</b>  | Anti-CCDC85C rabbit polyclonal antibody |
| <b>Applications</b> | ELISA, WB, IHC                          |
| <b>Immunogen</b>    | Synthetic peptide of human CCDC85C      |
| <b>Reactivity</b>   | Human, Mouse                            |
| <b>Content</b>      | 0.1 mg/ml                               |
| <b>Host species</b> | Rabbit                                  |
| <b>Ig class</b>     | Immunogen-specific rabbit IgG           |
| <b>Purification</b> | Antigen affinity purification           |

### Target information

|                  |                                   |
|------------------|-----------------------------------|
| <b>Symbol</b>    | CCDC85C                           |
| <b>Full name</b> | coiled-coil domain containing 85C |
| <b>Synonyms</b>  |                                   |
| <b>Swissprot</b> | A6NKD9                            |

### Target Background

CCDC85C may play an important role in cortical development, especially in the maintenance of radial glia.

订购热线: 4008-898-798


## Applications

### Immunohistochemistry

Predicted cell location: Cell junction

Positive control: Human liver cancer

Recommended dilution: 20-100



The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using ml163389(CCDC85C Antibody) at dilution 1/20, on the right is treated with synthetic peptide. (Original magnification:  $\times 200$ )

### Western blotting

Predicted band size: 45 kDa

Positive control: Human breast cancer tissue lysate

Recommended dilution: 200-1000

Gel: 8%SDS-PAGE


Lysate: 40  $\mu$ g

Lane: Human breast cancer tissue lysate

Primary antibody: ml163389(CCDC85C Antibody) at dilution 1/250

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 5 minutes



### ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: [mlbio\\_cn@yeah.net](mailto:mlbio_cn@yeah.net)

网址: [www.mlbio.cn](http://www.mlbio.cn)