

## Anti-NKX3-2 antibody

<b>Cat. No.</b>	ml163263
<b>Package</b>	25 µl/100 µl/200 µl
<b>Storage</b>	-20°C, pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol

### Product overview

<b>Description</b>	Anti-NKX3-2 rabbit polyclonal antibody
<b>Applications</b>	ELISA, IHC
<b>Immunogen</b>	Synthetic peptide of human NKX3-2
<b>Reactivity</b>	Human, Mouse
<b>Content</b>	0.2 mg/ml
<b>Host species</b>	Rabbit
<b>Ig class</b>	Immunogen-specific rabbit IgG
<b>Purification</b>	Antigen affinity purification

### Target information

<b>Symbol</b>	NKX3-2
<b>Full name</b>	NK3 homeobox 2
<b>Synonyms</b>	SMMD; BAPX1; NKX3B; NKX3.2
<b>Swissprot</b>	P78367

### Target Background

This gene encodes a member of the NK family of homeobox-containing proteins. The encoded protein may play a role in skeletal development.

订购热线: 4008-898-798


## Applications

### Immunohistochemistry


Predicted cell location: Nucleus

Positive control: Human thyroid cancer

Recommended dilution: 20-100




The image on the left is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using ml163263(NKX3-2 Antibody) at dilution 1/25, on the right is treated with synthetic peptide. (Original magnification:  $\times 200$ )




Predicted cell location: Nucleus

Positive control: Human esophagus cancer

Recommended dilution: 20-100



The image on the left is immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using ml163263(NKX3-2 Antibody) at dilution 1/25, on the right is treated with synthetic peptide. (Original magnification:  $\times 200$ )



### ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: [mlbio\\_cn@yeah.net](mailto:mlbio_cn@yeah.net)

网址: [www.mlbio.cn](http://www.mlbio.cn)