

订购热线: 4008-898-798

# Anti-FAM91A1 antibody

**Cat. No.** ml163882

Package 25  $\mu$ l/100  $\mu$ l/200  $\mu$ l

**Storage** -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

**Product overview** 

**Description** Anti-FAM91A1 rabbit polyclonal antibody

**Applications** ELISA, IHC

Immunogen Synthetic peptide of human FAM91A1

ReactivityHuman, MouseContent0.36 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

**Target information** 

Symbol FAM91A1

Full name family with sequence similarity 91 member A1

Synonyms

Swissprot Q658Y4

### **Target Background**

FAM91A1 (family with sequence similarity 91, member A1) is a 94 kDa phospho-protein that is not yet characterized.




订购热线: 4008-898-798


## **Applications**

#### **Immunohistochemistry**

Predicted cell location: Cytoplasm Positive control: Human thyroid cancer

Recommended dilution: 30-150







The image on the left is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using ml163882(FAM91A1 Antibody) at dilution 1/20, on the right is treated with synthetic peptide. (Original magnification: ×200)

Predicted cell location: Cytoplasm Positive control: Human cervical cancer

Recommended dilution: 30-150





The image on the left is immunohistochemistry of paraffin-embedded Human cervical cancer tissue using ml163882(FAM91A1 Antibody) at dilution 1/20, on the right is treated with synthetic peptide. (Original magnification: ×200)

#### **ELISA**

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio\_cn@yeah.net 网址: www.mlbio.cn