

订购热线: 4008-898-798

Anti-HMBOX1 antibody

Cat. No. ml164193

Package 25 μl/100 μl/200 μl

Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview

Description Anti-HMBOX1 rabbit polyclonal antibody

Applications ELISA, IHC

Immunogen Synthetic peptide of human HMBOX1

ReactivityHuman, MouseContent0.9 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

Target information

Symbol HMBOX1

Full namehomeobox containing 1SynonymsHOT1; PBHNF; HNF1LA

Swissprot Q6NT76

Target Background

Transcription factor. Isoform 1 acts as a transcriptional repressor. Isoform 4 has very low activity as a transcriptional repressor.



订购热线: 4008-898-798


Applications

Immunohistochemistry

Predicted cell location: Nucleus and Cytoplasm

Positive control: Human tonsil Recommended dilution: 30-150





The image on the left is immunohistochemistry of paraffin-embedded Human tonsil tissue using ml164193(HMBOX1 Antibody) at dilution 1/30, on the right is treated with synthetic peptide. (Original magnification: ×200)

Predicted cell location: Nucleus and Cytoplasm

Positive control: Human cervical cancer

Recommended dilution: 30-150





The image on the left is immunohistochemistry of paraffin-embedded Human cervical cancer tissue using ml164193(HMBOX1 Antibody) at dilution 1/30, on the right is treated with synthetic peptide. (Original magnification: ×200)

ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net 网址: www.mlbio.cn