

Anti-TBX6 antibody

 Cat. No.
 ml220640

 Package
 25 μl/100 μl/200 μl

 Storage
 -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview		
Description	Anti-TBX6 rabbit polyclonal an	tibody
Applications	ELISA, IHC	
Immunogen	Fusion protein of human TBX6	5
Reactivity	Human, Mouse, Rat	
Content	0.1 mg/ml	
Host species	Rabbit	
lg class	Immunogen-specific rabbit IgG	
Purification	Antigen affinity purification	444
Target information	(
Symbol	TBX6	
Full name	T-box 6	

SCDO5

O95947

70000

Target Background

Synonyms


Swissprot

This gene is a member of a phylogenetically conserved family of genes that share a common DNA-binding domain, the Tbox. T-box genes encode transcription factors involved in the regulation of developmental processes. Knockout studies in mice indicate that this gene is important for specification of paraxial mesoderm structures.



订购热线: 4008-898-798

Applications Immunohistochemistry Predicted cell location: Nucleus Positive control: Human colorectal cancer Recommended dilution: 20-100



The image is immunohistochemistry of paraffin-embedded Human colorectal cancer tissue using ml220640(TBX6 Antibody) at dilution 1/30. (Original magnification: ×200)

ELISA

Recommended dilution: 2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

- 邮箱:mlbio_cn@yeah.net 网址:www.mlbio.cn