

订购热线: 4008-898-798

# Anti-UGCG antibody

**Cat. No.** ml220685

Package 25  $\mu$ l/100  $\mu$ l/200  $\mu$ l

**Storage** -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

**Product overview** 

**Description** Anti-UGCG rabbit polyclonal antibody

Applications ELISA, IHC

Immunogen Fusion protein of human UGCG

Reactivity Human, Mouse, Rat

Content0.8 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

**Target information** 

Symbol UGCG

Full name UDP-glucose ceramide glucosyltransferase

Synonyms GCS; GLCT1
Swissprot Q16739

### **Target Background**

Glycosphingolipids (GSLs) are a group of membrane components that contain lipid and sugar moieties. They are present in essentially all animal cells and are believed to have important roles in various cellular processes. UDP-glucose ceramide glucosyltransferase catalyzes the first glycosylation step in glycosphingolipid biosynthesis. The product, glucosylceramide, is the core structure of more than 300 GSLs. UGCG is widely expressed and transcription is upregulated during keratinocyte differentiation.




订购热线: 4008-898-798


### **Applications**

## **Immunohistochemistry**

Predicted cell location: Cytoplasm Positive control: Human thyroid cancer

Recommended dilution: 50-200





The image on the left is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using ml220685(UGCG Antibody) at dilution 1/50, on the right is treated with fusion protein. (Original magnification: ×200)

#### **ELISA**

Recommended dilution: 2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio\_cn@yeah.net 网址: www.mlbio.cn