

订购热线: 4008-898-798

Anti-YY1AP1 antibody

Cat. No. ml223591

Package 25 μl/100 μl/200 μl

Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview

Description Anti-YY1AP1 rabbit polyclonal antibody

Applications ELISA, IHC

Immunogen Fusion protein of human YY1AP1

ReactivityHumanContent0.3 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

Target information

Symbol YY1AP1

Full name YY1 associated protein 1

Synonyms YAP; HCCA1; HCCA2; YY1AP

Swissprot Q9H869

Target Background

The encoded gene product presumably interacts with YY1 protein; however, its exact function is not known. Alternative splicing results in multiple transcript variants encoding different isoforms.




订购热线: 4008-898-798


Applications

Immunohistochemistry

Predicted cell location: Cytoplasm and Nucleus Positive control: Human cervical cancer

Recommended dilution: 20-100






The image on the left is immunohistochemistry of paraffin-embedded Human cervical cancer tissue using ml223591(YY1AP1 Antibody) at dilution 1/20, on the right is treated with fusion protein. (Original magnification: ×200)

Predicted cell location: Cytoplasm and Nucleus

Positive control: Human brain Recommended dilution: 20-100





The image on the left is immunohistochemistry of paraffin-embedded Human brain tissue using ml223591(YY1AP1 Antibody) at dilution 1/20, on the right is treated with fusion protein. (Original magnification: ×200)

ELISA

Recommended dilution: 2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net 网址: www.mlbio.cn