

Anti-THEMIS2 antibody

Cat. No.	ml223738
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-THEMIS2 rabbit polyclonal antibody
Applications	ELISA, IHC
Immunogen	Fusion protein of human THEMIS2
Reactivity	Human
Content	0.7 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	THEMIS2
Full name	thymocyte selection associated family member 2
Synonyms	ICB-1; C1orf38
Swissprot	Q5TEJ8

Target Background

May constitute a control point in macrophage inflammatory response, promoting LPS-induced TNF production.

订购热线: 4008-898-798


Applications

Immunohistochemistry

Predicted cell location: Cytoplasm and Cell membrane

Positive control: Human brain

Recommended dilution: 20-100



The image on the left is immunohistochemistry of paraffin-embedded Human brain tissue using ml223738(THEMIS2 Antibody) at dilution 1/20, on the right is treated with fusion protein. (Original magnification: $\times 200$)

ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn