

Anti-UBR2 antibody

Cat. No.	ml224155
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-UBR2 rabbit polyclonal antibody
Applications	ELISA, IHC
Immunogen	Fusion protein of human UBR2
Reactivity	Human, Mouse
Content	0.9 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	UBR2
Full name	ubiquitin protein ligase E3 component n-recognin 2
Synonyms	C6orf133; bA49A4.1; dJ242G1.1; dJ392M17.3
Swissprot	Q8I WV8

Target Background

This gene encodes an E3 ubiquitin ligase of the N-end rule proteolytic pathway that targets proteins with destabilizing N-terminal residues for polyubiquitylation and proteasome-mediated degradation. Alternative splicing results in multiple transcript variants.

订购热线: 4008-898-798


Applications

Immunohistochemistry

Predicted cell location: Cytoplasm and Nucleus

Positive control: Human lung cancer

Recommended dilution: 40-250



The image on the left is immunohistochemistry of paraffin-embedded Human lung cancer tissue using ml224155(UBR2 Antibody) at dilution 1/55, on the right is treated with fusion protein. (Original magnification: $\times 200$)

ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn