

订购热线: 4008-898-798

# Anti-USP35 antibody

**Cat. No.** ml224192

Package 25  $\mu$ l/100  $\mu$ l/200  $\mu$ l

**Storage** -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

**Product overview** 

**Description** Anti-USP35 rabbit polyclonal antibody

**Applications** ELISA, IHC

Immunogen Fusion protein of human USP35

ReactivityHumanContent0.4 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

**Target information** 

Symbol USP35

Full name ubiquitin specific peptidase 35

**Synonyms** 

Swissprot Q9P2H5

### **Target Background**

USP35 (ubiquitin specific peptidase 35) is a 1,017 amino acid protein that belongs to the peptidase C19 family. Expressed in pancreas, testis and skeletal muscle, USP35 functions to catalyze the conversion of a ubiquitin C-terminal thioester and water to free ubiquitin and a thiol, a reaction that may play a role in signaling events throughout the cell.




订购热线: 4008-898-798

### **Applications**

## **Immunohistochemistry**

Predicted cell location: Nucleus and Cytoplasm

Positive control: Human brain Recommended dilution: 20-100







The image on the left is immunohistochemistry of paraffin-embedded Human brain tissue using ml224192(USP35 Antibody) at dilution 1/30, on the right is treated with fusion protein. (Original magnification: ×200)

Predicted cell location: Nucleus and Cytoplasm

Positive control: Human prostate cancer



Recommended dilution: 20-100



The image on the left is immunohistochemistry of paraffin-embedded Human prostate cancer tissue using ml224192(USP35 Antibody) at dilution 1/30, on the right is treated with fusion protein. (Original magnification: ×200)

### **ELISA**

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio\_cn@yeah.net 网址: www.mlbio.cn