

Anti-ADAP2 antibody

Cat. No.	ml224508
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-ADAP2 rabbit polyclonal antibody
Applications	ELISA, IHC
Immunogen	Fusion protein of human ADAP2
Reactivity	Human, Mouse, Rat
Content	0.9 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	ADAP2
Full name	ArfGAP with dual PH domains 2
Synonyms	CENTA2; cent-b; HSA272195
Swissprot	Q9NPF8

Target Background

GTPase-activating protein for the ADP ribosylation factor family (Potential). Binds phosphatidylinositol 3,4,5-trisphosphate (PtdInsP₃) and inositol 1,3,4,5-tetrakisphosphate (InsP₄). Possesses a stoichiometry of two binding sites for InsP₄ with identical affinity.

订购热线: 4008-898-798



Applications

Immunohistochemistry

Predicted cell location: Cytoplasm and Cell membrane

Positive control: Human liver cancer

Recommended dilution: 25-100



The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using ml224508(ADAP2 Antibody) at dilution 1/20, on the right is treated with fusion protein. (Original magnification: $\times 200$)

ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn