

Anti-CFHR2 antibody

 Cat. No.
 ml224525

 Package
 25 μl/100 μl/200 μl

 Storage
 -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview	
Description	Anti-CFHR2 rabbit polyclonal antibody
Applications	ELISA, IHC
Immunogen	Fusion protein of human CFHR2
Reactivity	Human
Content	0.5 mg/ml
Host species	Rabbit
lg class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification
Target information	
Symbol	CFHR2

noc

Symbol Full name Synonyms Swissprot CFHR2 complement factor H related 2 FHR2; HFL3; CFHL2 P36980

Target Background


This gene belongs to a family of complement factor H-related genes (CFHR), which are clustered together with complement factor H gene on chromosome 1, and are involved in regulation of complement. Mutations in CFHR genes have been associated with dense deposit disease and atypical haemolytic-uraemic syndrome. Alternatively spliced transcript variants have been found for this gene.



订购热线: 4008-898-798

Applications Immunohistochemistry

Predicted cell location: Cytoplasm or Nucleus Positive control: Human gastric cancer Recommended dilution: 25-100




The image on the left is immunohistochemistry of paraffin-embedded Human gastric cancer tissue using ml224525(CFHR2 Antibody) at dilution 1/25, on the right is treated with fusion protein. (Original magnification: ×200)

ELISA

Recommended dilution: 5000-10000

Predicted cell location: Cytoplasm or Nucleus Positive control: Human colorectal cancer Recommended dilution: 25-100





The image on the left is immunohistochemistry of paraffin-embedded Human colorectal cancer tissue using ml224525(CFHR2 Antibody) at dilution 1/25, on the right is treated with fusion protein. (Original magnification: ×200)

- 联系电话: 4008-898-798, 021-61725725
- 联系QQ: 2881505695,2881505696
- 邮箱: mlbio_cn@yeah.net 网址: www.mlbio.cn