

订购热线: 4008-898-798

Anti-CHAC2 antibody

Cat. No.	ml224530
Package	25 μl/100 μl/200 μl
Storage	-20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview	
Description	Anti-CHAC2 rabbit polyclonal antibody
Applications	ELISA, IHC
Immunogen	Full length fusion protein
Reactivity	Human, Mouse, Rat
Content	0.5 mg/ml
Host species	Rabbit
lg class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification
Target information	
Symbol	CHAC2
Full name	ChaC cation transport regulator homolog 2
Synonyms	
Swissprot	Q8WUX2

300del

Target Background

The protein encoded by this gene is a gamma-glutamyl cyclotransferase that catalyzes the conversion of glutathione to 5oxoproline and cysteinylglycine. It is thought that this gene is upregulated in response to endoplasmic reticulum stress and that the glutathione depletion enhances apoptosis.


duce



订购热线: 4008-898-798

Applications Immunohistochemistry

Predicted cell location: Cytoplasm and Cell membrane Positive control: Human liver cancer Recommended dilution: 25-100



The image is immunohistochemistry of paraffin-embedded Human liver cancer tissue using ml224530(CHAC2 Antibody) at dilution 1/20. (Original magnification: ×200)

ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱:mlbio_cn@yeah.net 网址:www.mlbio.cn