

Anti-CLUL1 antibody

Cat. No.	ml224592
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-CLUL1 rabbit polyclonal antibody
Applications	ELISA, IHC
Immunogen	Fusion protein of human CLUL1
Reactivity	Human
Content	0.5 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	CLUL1
Full name	clusterin like 1
Synonyms	RA337M
Swissprot	Q15846

Target Background

CLUL1 is a retina-specific clusterin-like protein which shows nearly 25% identity to clusterin at the protein level. It is expressed specifically in cone photoreceptor cells.

订购热线: 4008-898-798


Applications

Immunohistochemistry

Predicted cell location: Secreted

Positive control: Human liver cancer

Recommended dilution: 20-100



The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using ml224592(CLUL1 Antibody) at dilution 1/20, on the right is treated with fusion protein. (Original magnification: $\times 200$)

ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn