

Anti-CHKB-CPT1B antibody

 Cat. No.
 ml224669

 Package
 25 μl/100 μl/200 μl

 Storage
 -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview	
Description	Anti-CHKB-CPT1B rabbit polyclonal antibody
Applications	ELISA, IHC
Immunogen	Fusion protein of human CHKB-CPT1B
Reactivity	Human, Mouse, Rat
Content	0.6 mg/ml
Host species	Rabbit
lg class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification
Target information	marks.
Symbol	СНКВ-СРТ1В

CHKB-CPT1B CHKB-CPT1B readthrough (NMD candidate) CHKL-CPT1B Q92523

, cê

Target Background

Full name

Synonyms


Swissprot

The genes CHKB and CPT1B are adjacent on chromosome 22 and read-through transcripts are expressed that include exons from both loci. The read-through transcripts are candidates for nonsense-mediated mRNA decay (NMD) and are unlikely to express proteins.




Applications Immunohistochemistry


Predicted cell location: Cytoplasm Positive control: Human liver cancer Recommended dilution: 25-100



The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using ml224669(CHKB-CPT1B Antibody) at dilution 1/25, on the right is treated with fusion protein. (Original magnification: ×200)

Predicted cell location: Cytoplasm Positive control: Human gastric cancer Recommended dilution: 25-100





The image on the left is immunohistochemistry of paraffin-embedded Human gastric cancer tissue using ml224669(CHKB-CPT1B Antibody) at dilution 1/25, on the right is treated with fusion protein. (Original magnification: ×200)

Recommended dilution: 5000-10000

ELISA

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695,2881505696、

邮箱: mlbio_cn@yeah.net 网址: www.mlbio.cn