

## Anti-SGCA antibody

 Cat. No.
 ml224784

 Package
 25 μl/100 μl/200 μl

 Storage
 -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview	
Description	Anti-SGCA rabbit polyclonal antibody
Applications	ELISA, IHC
Immunogen	Fusion protein of human SGCA
Reactivity	Human, Mouse
Content	1.2 mg/ml
Host species	Rabbit
lg class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification
Target information	
Symbol	SGCA
Full name	sarcoglycan alpha
Synonyms	ADL; DAG2; 50DAG; DMDA2; LGMD2D; SCARMD1; adhalin
Swissprot	Q16586

## **Target Background**

1000

This gene encodes a component of the dystrophin-glycoprotein complex (DGC), which is critical to the stability of muscle fiber membranes and to the linking of the actin cytoskeleton to the extracellular matrix. Its expression is thought to be restricted to striated muscle. Mutations in this gene result in type 2D autosomal recessive limb-girdle muscular dystrophy. Multiple transcript variants encoding different isoforms have been found for this gene.


ducer



订购热线: 4008-898-798

## Applications Immunohistochemistry

Predicted cell location: Cytoplasm or Nucleus Positive control: Human liver cancer Recommended dilution: 25-100





The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using ml224784(SGCA Antibody) at dilution 1/35, on the right is treated with fusion protein. (Original magnification: ×200)

## ELISA

Recommended dilution: 5000-10000

Predicted cell location: Cytoplasm or Nucleus Positive control: Human lung cancer Recommended dilution: 25-100





The image on the left is immunohistochemistry of paraffin-embedded Human lung cancer tissue using ml224784(SGCA Antibody) at dilution 1/35, on the right is treated with fusion protein. (Original magnification: ×200)

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695,2881505696

邮箱: mlbio\_cn@yeah.net 网址: www.mlbio.cn