

Anti-PDRG1 antibody

Cat. No.	ml225237
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-PDRG1 rabbit polyclonal antibody
Applications	ELISA, IHC
Immunogen	Fusion protein of human PDRG1
Reactivity	Human, Mouse, Rat
Content	1.32 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	PDRG1
Full name	p53 and DNA damage regulated 1
Synonyms	PDRG; C20orf126
Swissprot	Q9NUG6

Target Background

May play a role in chaperone-mediated protein folding.

订购热线: 4008-898-798


Applications

Immunohistochemistry

Predicted cell location: Cytoplasm

Positive control: Human liver cancer

Recommended dilution: 50-300



The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using ml225237(PDRG1 Antibody) at dilution 1/75, on the right is treated with fusion protein. (Original magnification: ×200)

ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn