

Anti-THAP12 antibody

Cat. No.	ml225250
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-THAP12 rabbit polyclonal antibody
Applications	ELISA, IHC
Immunogen	Fusion protein of human THAP12
Reactivity	Human, Mouse
Content	0.78 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	THAP12
Full name	THAP domain containing 12
Synonyms	DAP4; THAP0; PRKRIR; P52rIPK
Swissprot	O43422

Target Background

Upstream regulator of interferon-induced serine/threonine protein kinase R (PKR). May block the PKR-inhibitory function of DNAJC3, resulting in restoration of kinase activity and suppression of cell growth.

订购热线: 4008-898-798


Applications

Immunohistochemistry

Predicted cell location: Nucleus

Positive control: Human ovarian cancer

Recommended dilution: 40-200



The image is immunohistochemistry of paraffin-embedded Human ovarian cancer tissue using ml225250(THAP12 Antibody) at dilution 1/45. (Original magnification: ×200)

ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn