

Anti-GSG1L antibody

Cat. No.	ml225611
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-GSG1L rabbit polyclonal antibody
Applications	ELISA, IHC
Immunogen	Fusion protein of human GSG1L
Reactivity	Human, Mouse, Rat
Content	0.9 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	GSG1L
Full name	GSG1 like
Synonyms	PRO19651
Swissprot	Q6UXU4

Target Background

As a component of the inner core of AMPAR complex, modifies AMPA receptor (AMPA) gating.



Applications

Immunohistochemistry

Predicted cell location: Cytoplasm and Nucleus

Positive control: Human breast cancer

Recommended dilution: 50-300





The image on the left is immunohistochemistry of paraffin-embedded Human breast cancer tissue using ml225611(GSG1L Antibody) at dilution 1/50, on the right is treated with fusion protein. (Original magnification: $\times 200$)

Predicted cell location: Cytoplasm and Nucleus

Positive control: Human ovarian cancer

Recommended dilution: 50-300



The image on the left is immunohistochemistry of paraffin-embedded Human ovarian cancer tissue using ml225611(GSG1L Antibody) at dilution 1/50, on the right is treated with fusion protein. (Original magnification: $\times 200$)

ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn