

Anti-PPM1B antibody

Cat. No. ml225633 Package 25 µl/100 µl/200 µl -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol Storage

Product overview	
Description	Anti-PPM1B rabbit polyclonal antibody
Applications	ELISA, IHC
Immunogen	Fusion protein of human PPM1B
Reactivity	Human
Content	0.72 mg/ml
Host species	Rabbit
lg class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification
Target information	
Symbol	PPM1B
Full name	protein phosphatase, Mg2+/Mn2+ dependent 1B

Symbol	PPM1B
Full name	protein phosphatase, Mg2+/Mn2+ dependent 1B
Synonyms	PP2CB; PP2CBETA; PP2C-beta; PPC2BETAX; PP2C-beta-X
Swissprot	O75688
Target Background	du


Target Background

The protein encoded by this gene is a member of the PP2C family of Ser/Thr protein phosphatases. PP2C family members are known to be negative regulators of cell stress response pathways. This phosphatase has been shown to dephosphorylate cyclin-dependent kinases (CDKs), and thus may be involved in cell cycle control. Overexpression of this phosphatase is reported to cause cell-growth arrest or cell death. Alternative splicing results in multiple transcript variants encoding different isoforms. Additional transcript variants have been described, but currently do not represent full-length sequences.



订购热线: 4008-898-798

Applications Immunohistochemistry Predicted cell location: Cytoplasm Positive control: Human tonsil Recommended dilution: 50-300



The image on the left is immunohistochemistry of paraffin-embedded Human tonsil tissue using ml225633(PPM1B Antibody) at dilution 1/50, on the right is treated with fusion protein. (Original magnification: ×200)

ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695,2881505696

邮箱: mlbio_cn@yeah.net 网址: www.mlbio.cn