

Anti-DCAF10 antibody

Cat. No.	ml225666
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-DCAF10 rabbit polyclonal antibody
Applications	ELISA, IHC
Immunogen	Fusion protein of human DCAF10
Reactivity	Human, Mouse
Content	2.28 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	DCAF10
Full name	DDB1 and CUL4 associated factor 10
Synonyms	WDR32
Swissprot	Q5QP82

Target Background

May function as a substrate receptor for CUL4-DDB1 E3 ubiquitin-protein ligase complex.



Applications

Immunohistochemistry

Predicted cell location: Nucleus

Positive control: Human thyroid cancer

Recommended dilution: 50-300





The image on the left is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using ml225666(DCAF10 Antibody) at dilution 1/130, on the right is treated with fusion protein. (Original magnification: $\times 200$)

Predicted cell location: Nucleus

Positive control: Human esophagus cancer

Recommended dilution: 50-300



The image on the left is immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using ml225666(DCAF10 Antibody) at dilution 1/130, on the right is treated with fusion protein. (Original magnification: $\times 200$)

ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn