

Anti-RTKN2 antibody

Cat. No.	ml225757
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-RTKN2 rabbit polyclonal antibody
Applications	ELISA, IHC
Immunogen	Fusion protein of human RTKN2
Reactivity	Human
Content	1.32 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	RTKN2
Full name	rhotekin 2
Synonyms	PLEKHK1; bA531F24.1
Swissprot	Q8IZC4

Target Background

May play an important role in lymphopoiesis.

订购热线: 4008-898-798



Applications

Immunohistochemistry

Predicted cell location: Cytoplasm

Positive control: Human colorectal cancer

Recommended dilution: 50-300





The image on the left is immunohistochemistry of paraffin-embedded Human colorectal cancer tissue using ml225757(RTKN2 Antibody) at dilution 1/75, on the right is treated with fusion protein. (Original magnification: $\times 200$)

Predicted cell location: Cytoplasm

Positive control: Human liver cancer

Recommended dilution: 50-300



The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using ml225757(RTKN2 Antibody) at dilution 1/75, on the right is treated with fusion protein. (Original magnification: $\times 200$)

ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn