

Anti-IQCD antibody

Cat. No.	ml225770
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-IQCD rabbit polyclonal antibody
Applications	ELISA, IHC
Immunogen	Fusion protein of human IQCD
Reactivity	Human
Content	1.62 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	IQCD
Full name	IQ motif containing D
Synonyms	DRC10; CFAP84; 4933433C09Rik
Swissprot	Q96DY2

Target Background

IQCD contains 1 IQ domain. The function of the IQCD protein is not known.

订购热线: 4008-898-798


Applications

Immunohistochemistry

Predicted cell location: Cytoplasm

Positive control: Human colorectal cancer

Recommended dilution: 50-300



The image on the left is immunohistochemistry of paraffin-embedded Human colorectal cancer tissue using ml225770(IQCD Antibody) at dilution 1/90, on the right is treated with fusion protein. (Original magnification: $\times 200$)

ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn