

Anti-TP53BP2 antibody

 Cat. No.
 ml260110

 Package
 25 μl/100 μl/200 μl

 Storage
 -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview	
Description	Anti-TP53BP2 rabbit polyclonal antibody
Applications	ELISA, IHC
Immunogen	Synthetic peptide of human TP53BP2
Reactivity	Human, Mouse
Content	0.1 mg/ml
Host species	Rabbit
lg class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification
Target information	
Symbol	TP53BP2

TP53BP2 tumor protein p53 binding protein, 2 BBP, 53BP2, ASPP2, P53BP2, PPP1R13A Q13625

Target Background

Full name

Synonyms Swissprot


This gene encodes a member of the ASPP (apoptosis-stimulating protein of p53) family of p53 interacting proteins. The protein contains four ankyrin repeats and an SH3 domain involved in protein-protein interactions. It is localized to the perinuclear region of the cytoplasm, and regulates apoptosis and cell growth through interactions with other regulatory molecules including members of the p53 family. Multiple transcript variants encoding different isoforms have been found for this gene.



订购热线: 4008-898-798

Applications Immunohistochemistry

Predicted cell location: Cytoplasm Positive control: Human gasrtic cancer Recommended dilution: 5-20




The image on the left is immunohistochemistry of paraffin-embedded Human gasrtic cancer tissue using ml260110(TP53BP2 Antibody) at dilution 1/5, on the right is treated with synthetic peptide. (Original magnification: ×200)

ELISA

Recommended dilution: 1000-2000

Predicted cell location: Cytoplasm Positive control: Human colon cancer Recommended dilution: 5-20





The image on the left is immunohistochemistry of paraffin-embedded Human colon cancer tissue using ml260110(TP53BP2 Antibody) at dilution 1/5, on the right is treated with synthetic peptide. (Original magnification: ×200)

- 联系电话: 4008-898-798, 021-61725725
- 联系QQ: 2881505695,2881505696、
- 邮箱: mlbio_cn@yeah.net 网址: www.mlbio.cn