

订购热线: 4008-898-798

Anti-FGFR3 antibody

Cat. No. ml260329

Package 25 μ l/100 μ l/200 μ l

Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview

Description Anti-FGFR3 rabbit polyclonal antibody

Applications ELISA, IHC

Immunogen Synthetic peptide of human FGFR3

Reactivity Human, Mouse
Content 0.3 mg/ml
Host species Rabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

Target information

Symbol FGFR3

Full name Fibroblast growth factor receptor 3

Synonyms ACH, CEK2, JTK4, CD333, HSFGFR3EX

Swissprot P22607

Target Background

This gene encodes a member of the fibroblast growth factor receptor (FGFR) family, with its amino acid sequence being highly conserved between members and among divergent species. FGFR family members differ from one another in their ligand affinities and tissue distribution. A full-length representative protein would consist of an extracellular region, composed of three immunoglobulin-like domains, a single hydrophobic membrane-spanning segment and a cytoplasmic tyrosine kinase domain. The extracellular portion of the protein interacts with fibroblast growth factors, setting in motion a cascade of downstream signals, ultimately influencing mitogenesis and differentiation. This particular family member binds acidic and basic fibroblast growth hormone and plays a role in bone development and maintenance. Mutations in this gene lead to craniosynostosis and multiple types of skeletal dysplasia. Three alternatively spliced transcript variants that encode different protein isoforms have been described.




订购热线: 4008-898-798

Applications

Immunohistochemistry


Predicted cell location: Cytoplasm Positive control: Human liver cancer Recommended dilution: 35-150






The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using ml260329(FGFR3 Antibody) at dilution 1/35, on the right is treated with synthetic peptide. (Original magnification: ×200)

Predicted cell location: Cytoplasm Positive control: Human breast cancer Recommended dilution: 35-150





The image on the left is immunohistochemistry of paraffin-embedded Human breast cancer tissue using ml260329(FGFR3 Antibody) at dilution 1/35, on the right is treated with synthetic peptide. (Original magnification: ×200)

ELISA

Recommended dilution: 1000-2000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net 网址: www.mlbio.cn