

订购热线: 4008-898-798

Anti-PTPRT antibody

Cat. No. ml261100

Package 25 μl/100 μl/200 μl

Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview

Description Anti-PTPRT rabbit polyclonal antibody

Applications ELISA, IHC

Immunogen Synthetic peptide of human PTPRT

ReactivityHuman, MouseContent0.5 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

Target information

Symbol PTPRT

Full name Protein tyrosine phosphatase, receptor type, T

Synonyms RPTPrho
Swissprot 014522

Target Background

The protein encoded by this gene is a member of the protein tyrosine phosphatase (PTP) family. PTPs are known to be signaling molecules that regulate a variety of cellular processes including cell growth, differentiation, mitotic cycle, and oncogenic transformation. This PTP possesses an extracellular region, a single transmembrane region, and two tandem intracellular catalytic domains, and thus represents a receptor-type PTP. The extracellular region contains a meprin-A5 antigen-PTP (MAM) domain, Ig-like and fibronectin type III-like repeats. The protein domain structure and the expression pattern of the mouse counterpart of this PTP suggest its roles in both signal transduction and cellular adhesion in the central nervous system. Two alternatively spliced transcript variants of this gene, which encode distinct proteins, have been reported.




订购热线: 4008-898-798


Applications

Immunohistochemistry

Predicted cell location: Cytoplasm Positive control: Human ovarian cancer

Recommended dilution: 25-100





The image on the left is immunohistochemistry of paraffin-embedded Human ovarian cancer tissue using ml261100(PTPRT Antibody) at dilution 1/30, on the right is treated with synthetic peptide. (Original magnification: ×200)

ELISA

Recommended dilution: 1000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net 网址: www.mlbio.cn