

订购热线: 4008-898-798

# Anti-NPAP1 antibody

**Cat. No.** ml261393

**Package** 25 μl/100 μl/200 μl

**Storage** -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

**Product overview** 

**Description** Anti-NPAP1 rabbit polyclonal antibody

Applications ELISA, IHC

Immunogen Synthetic peptide of human NPAP1

ReactivityHumanContent1.6 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

**Target information** 

Symbol NPAP1

Full name nuclear pore associated protein 1

SynonymsC15orf2SwissprotQ9NZP6

### **Target Background**

This gene is located in the Prader-Willi syndrome region on chromosome 15. This gene is biallelically expressed in adult testis and brain but is paternally imprinted in fetal brain. Defects in this gene may be associated with Prader-Willi syndrome.




订购热线: 4008-898-798


# **Applications**

## **Immunohistochemistry**

Predicted cell location: Cytoplasm Positive control: Human thyroid cancer

Recommended dilution: 50-200





The image on the left is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using ml261393(NPAP1 Antibody) at dilution 1/70, on the right is treated with synthetic peptide. (Original magnification: ×200)

#### **ELISA**

Recommended dilution: 1000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio\_cn@yeah.net 网址: www.mlbio.cn