

订购热线: 4008-898-798

# Anti-HP antibody

**Cat. No.** ml261755

Package 25  $\mu$ l/100  $\mu$ l/200  $\mu$ l

**Storage** -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

**Product overview** 

**Description** Anti-HP rabbit polyclonal antibody

**Applications** ELISA, IHC

Immunogen Synthetic peptide of human HP

ReactivityHumanContent0.5 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

**Target information** 

Symbol HP

Full name haptoglobin

Synonyms BP; HPA1S; HP2ALPHA2

Swissprot P00738

### **Target Background**

This gene encodes a preproprotein, which is processed to yield both alpha and beta chains, which subsequently combine as a tetramer to produce haptoglobin. Haptoglobin functions to bind free plasma hemoglobin, which allows degradative enzymes to gain access to the hemoglobin, while at the same time preventing loss of iron through the kidneys and protecting the kidneys from damage by hemoglobin. Mutations in this gene and/or its regulatory regions cause ahaptoglobinemia or hypohaptoglobinemia.





订购热线: 4008-898-798

### **Applications**

## **Immunohistochemistry**


Predicted cell location: Cytoplasm Positive control: Human gasrtic cancer Recommended dilution: 100-300






The image on the left is immunohistochemistry of paraffin-embedded Human gasrtic cancer tissue using ml261755(HP Antibody) at dilution 1/40, on the right is treated with synthetic peptide. (Original magnification: ×200)

Predicted cell location: Cytoplasm Positive control: Human thyroid cancer Recommended dilution: 100-300





The image on the left is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using ml261755(HP Antibody) at dilution 1/40, on the right is treated with synthetic peptide. (Original magnification: ×200)

#### **ELISA**

Recommended dilution: 2000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio\_cn@yeah.net 网址: www.mlbio.cn