

订购热线: 4008-898-798

Anti-SCD antibody

Cat. No. ml262163

Package 25 μl/100 μl/200 μl

Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview

Description Anti-SCD rabbit polyclonal antibody

Applications ELISA, IHC

Immunogen Synthetic peptide of human SCD

ReactivityHumanContent0.5 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

Target information

Symbol SCD

Full name stearoyl-CoA desaturase (delta-9-desaturase)

Synonyms SCD1; FADS5; SCDOS; MSTP008

Swissprot 000767

Target Background

This gene encodes an enzyme involved in fatty acid biosynthesis, primarily the synthesis of oleic acid. Transcripts of approximately 3.9 and 5.2 kb, differing only by alternative polyadenlyation signals, have been detected. A gene encoding a similar enzyme is located on chromosome 4 and a pseudogene of this gene is located on chromosome 17.





订购热线: 4008-898-798

Applications

Immunohistochemistry


Predicted cell location: Cytoplasm Positive control: Human lung cancer Recommended dilution: 25-100






The image on the left is immunohistochemistry of paraffin-embedded Human lung cancer tissue using ml262163(SCD Antibody) at dilution 1/25, on the right is treated with synthetic peptide. (Original magnification: ×200)

Predicted cell location: Cytoplasm Positive control: Human brain Recommended dilution: 25-100





The image on the left is immunohistochemistry of paraffin-embedded Human brain tissue using ml262163(SCD Antibody) at dilution 1/25, on the right is treated with synthetic peptide. (Original magnification: ×200)

ELISA

Recommended dilution: 1000-2000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net 网址: www.mlbio.cn