

订购热线: 4008-898-798

# Anti-SEMA7A antibody

**Cat. No.** ml262182

**Package** 25 μl/100 μl/200 μl

**Storage** -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

**Product overview** 

**Description** Anti-SEMA7A rabbit polyclonal antibody

Applications ELISA, IHC

Immunogen Synthetic peptide of human SEMA7A

ReactivityHuman, MouseContent0.6 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

**Target information** 

Symbol SEMA7A

Full name semaphorin 7A, GPI membrane anchor (John Milton Hagen blood group)

Synonyms JMH; CD108; SEMAL; CDw108; SEMAK1; H-Sema-L; H-SEMA-K1

Swissprot 075326

## **Target Background**

The protein encoded by this gene binds to cell surfaces through a glycosylphosphatidylinositol (GPI) linkage. The encoded glycoprotein is found on activated lymphocytes and erythrocytes. This protein may be involved in immunomodulatory and neuronal processes. Defects in this gene can result in loss of bone mineral density (BMD). Three transcript variants encoding different isoforms have been found for this gene.




订购热线: 4008-898-798


### **Applications**

## **Immunohistochemistry**

Predicted cell location: Cytoplasm or Cell membrane

Positive control: Human lung cancer Recommended dilution: 25-100





The image on the left is immunohistochemistry of paraffin-embedded Human lung cancer tissue using ml262182(SEMA7A Antibody) at dilution 1/35, on the right is treated with synthetic peptide. (Original magnification: ×200)

#### **ELISA**

Recommended dilution: 1000-2000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio\_cn@yeah.net 网址: www.mlbio.cn