

订购热线: 4008-898-798

Anti-TCEB3B antibody

Cat. No. ml262345

Package 25 μ l/100 μ l/200 μ l

Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview

Description Anti-TCEB3B rabbit polyclonal antibody

Applications ELISA, IHC

Immunogen Synthetic peptide of human TCEB3B

ReactivityHumanContent0.3 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

Target information

Symbol TCEB3B

Full name transcription elongation factor B polypeptide 3B (elongin A2)

Synonyms ELOA2; HsT832; TCEB3L

Swissprot Q8IYF1

Target Background

This gene encodes the transcriptionally active subunit of the SIII (or elongin) transcription elongation factor complex, which also includes two regulatory subunits, elongins B and C. This complex acts to increase the rate of RNA chain elongation by RNA polymerase II by suppressing transient pausing of the polymerase at many sites along the DNA template. Whereas a related protein with similar function, elongin A, is ubiquitously expressed, the encoded protein is specifically expressed in the testis, suggesting it may have a role in spermatogenesis.




订购热线: 4008-898-798


Applications

Immunohistochemistry

Predicted cell location: Cytoplasm and Nucleus Positive control: Human esophagus cancer

Recommended dilution: 30-150







The image on the left is immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using ml262345(TCEB3B Antibody) at dilution 1/25, on the right is treated with synthetic peptide. (Original magnification: ×200)

Predicted cell location: Cytoplasm and Nucleus

Positive control: Human lung cancer Recommended dilution: 30-150





The image on the left is immunohistochemistry of paraffin-embedded Human lung cancer tissue using ml262345(TCEB3B Antibody) at dilution 1/25, on the right is treated with synthetic peptide. (Original magnification: ×200)

ELISA

Recommended dilution: 2000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net 网址: www.mlbio.cn