

订购热线: 4008-898-798

# Anti-TGM4 antibody

**Cat. No.** ml262352

**Package** 25 μl/100 μl/200 μl

**Storage** -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

**Product overview** 

**Description** Anti-TGM4 rabbit polyclonal antibody

Applications ELISA, IHC

Immunogen Synthetic peptide of human TGM4

ReactivityHumanContent0.9 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

**Target information** 

Symbol TGM4

Full name transglutaminase 4

SynonymsTGP; hTGPSwissprotP49221

## **Target Background**

Associated with the mammalian reproductive process. Catalyzes the cross-linking of proteins and the conjugation of polyamines to specific proteins in the seminal tract.




订购热线: 4008-898-798


#### **Applications**

### **Immunohistochemistry**

Predicted cell location: Cytoplasm Positive control: Human prostate cancer

Recommended dilution: 10-50





The image on the left is immunohistochemistry of paraffin-embedded Human prostate cancer tissue using ml262352(TGM4 Antibody) at dilution 1/20, on the right is treated with synthetic peptide. (Original magnification: ×200)

#### **ELISA**

Recommended dilution: 1000-2000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio\_cn@yeah.net 网址: www.mlbio.cn