

Anti-HAUS7 antibody

Cat. No.	ml262445
Package	25 μl/100 μl/200 μl
Storage	-20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview	
Description	Anti-HAUS7 rabbit polyclonal antibody
Applications	ELISA, IHC
Immunogen	Synthetic peptide of human HAUS7
Reactivity	Human
Content	0.6 mg/ml
Host species	Rabbit
lg class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification
Target information	
Symbol	HAUS7
Full name	HAUS augmin-like complex, subunit 7

30000

HAUS7 HAUS augmin-like complex, subunit 7 UIP1; UCHL5IP Q99871

Target Background

Synonyms

Swissprot


This gene encodes a subunit of the augmin complex, which regulates centrosome and mitotic spindle integrity, and is necessary for the completion of cytokinesis. The encoded protein was identified by interaction with ubiquitin C-terminal hydrolase 37. Alternative splicing results in multiple transcript variants.



订购热线: 4008-898-798

Applications Immunohistochemistry

Predicted cell location: Cytoplasm Positive control: Human liver cancer Recommended dilution: 25-100



The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using ml262445(HAUS7 Antibody) at dilution 1/30, on the right is treated with synthetic peptide. (Original magnification: ×200)

ELISA

Recommended dilution: 2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695,2881505696

邮箱: mlbio_cn@yeah.net 网址: www.mlbio.cn