

Anti-CHRM5 antibody

 Cat. No.
 ml262682

 Package
 25 μl/100 μl/200 μl

 Storage
 -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview	
Description	Anti-CHRM5 rabbit polyclonal antibody
Applications	ELISA, IHC
Immunogen	Synthetic peptide of human CHRM5
Reactivity	Human, Mouse, Rat
Content	0.7 mg/ml
Host species	Rabbit
lg class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification
Target information	
Symbol	CHRM5
Full name	cholinergic receptor, muscarinic 5
Synonyms	HM5

P08912

Target Background

Swissprot


The muscarinic cholinergic receptors belong to a larger family of G protein-coupled receptors. The functional diversity of these receptors is defined by the binding of acetylcholine and includes cellular responses such as adenylate cyclase inhibition, phosphoinositide degeneration, and potassium channel mediation. Muscarinic receptors influence many effects of acetylcholine in the central and peripheral nervous system. The clinical implications of this receptor are unknown; however, stimulation of this receptor is known to increase cyclic AMP levels.



订购热线: 4008-898-798

Applications Immunohistochemistry

Predicted cell location: Cell membrane Positive control: Human liver cancer Recommended dilution: 25-100





The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using ml262682(CHRM5 Antibody) at dilution 1/25, on the right is treated with synthetic peptide. (Original magnification: ×200)

ELISA

Recommended dilution: 2000-5000

Predicted cell location: Cell membrane Positive control: Human breast cancer Recommended dilution: 25-100





The image on the left is immunohistochemistry of paraffin-embedded Human breast cancer tissue using ml262682(CHRM5 Antibody) at dilution 1/25, on the right is treated with synthetic peptide. (Original magnification: ×200)

- 联系电话: 4008-898-798, 021-61725725
- 联系QQ: 2881505695,2881505696、
- 邮箱: mlbio_cn@yeah.net 网址: www.mlbio.cn