

订购热线: 4008-898-798

Anti-GJB5 antibody

Cat. No. ml262692

Package 25 μl/100 μl/200 μl

Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview

Description Anti-GJB5 rabbit polyclonal antibody

Applications ELISA, IHC

Immunogen Synthetic peptide of human GJB5

Reactivity Human, Mouse, Rat

Content 1.4 mg/ml Host species Rabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

Target information

Symbol GJB5

Full name gap junction protein, beta 5, 31.1kDa

Synonyms CX31.1 Swissprot O95377

Target Background

This gene encodes a member of the beta-type (group I) connexin family. The encoded protein is a gap junction protein involved in intercellular communication related to epidermal differentiation and environmental sensing. This gene has been linked to non-small cell lung cancer.





订购热线: 4008-898-798

Applications

Immunohistochemistry

Predicted cell location: Cytoplasm Positive control: Human brain Recommended dilution: 30-150





The image on the left is immunohistochemistry of paraffin-embedded Human brain tissue using ml262692(GJB5 Antibody) at dilution 1/45, on the right is treated with synthetic peptide. (Original magnification: ×200)

ELISA

Recommended dilution: 2000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn