

Anti-COX17 antibody

Cat. No. ml262710 Package 25 µl/100 µl/200 µl -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol Storage

Product overview	
Description	Anti-COX17 rabbit polyclonal antibody
Applications	ELISA, IHC
Immunogen	Synthetic peptide of human COX17
Reactivity	Human, Mouse
Content	1.2 mg/ml
Host species	Rabbit
lg class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification
Target information	
Symbol	COX17
Full name	COX17 cytochrome c oxidase copper chaperone
Synonyms	
Swissprot	Q14061
	JUCC.
Target Background	

Target Background

Cytochrome c oxidase (COX), the terminal component of the mitochondrial respiratory chain, catalyzes the electron transfer from reduced cytochrome c to oxygen. This component is a heteromeric complex consisting of 3 catalytic subunits encoded by mitochondrial genes and multiple structural subunits encoded by nuclear genes. The mitochondrially-encoded subunits function in electron transfer, and the nuclear-encoded subunits may function in the regulation and assembly of the complex. This nuclear gene encodes a protein which is not a structural subunit, but may be involved in the recruitment of copper to mitochondria for incorporation into the COX apoenzyme. This protein shares 92% amino acid sequence identity with mouse and rat Cox17 proteins. This gene is no longer considered to be a candidate gene for COX deficiency. A pseudogene COX17P has been found on chromosome 13.



订购热线: 4008-898-798

Applications Immunohistochemistry Predicted cell location: Cytoplasm Positive control: Human brain Recommended dilution: 30-150



The image on the left is immunohistochemistry of paraffin-embedded Human brain tissue using ml262710(COX17 Antibody) at dilution 1/45, on the right is treated with synthetic peptide. (Original magnification: ×200)

ELISA

Recommended dilution: 2000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695,2881505696

邮箱: mlbio_cn@yeah.net 网址: www.mlbio.cn