

订购热线: 4008-898-798

Anti-DVL3 antibody

Cat. No. ml262731

Package 25 μl/100 μl/200 μl

Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview

Description Anti-DVL3 rabbit polyclonal antibody

Applications ELISA, IHC

Immunogen Synthetic peptide of human DVL3

ReactivityHuman, MouseContent0.5 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

Target information

Symbol DVL3

Full name dishevelled segment polarity protein 3

Synonyms

Swissprot Q92997

Target Background

This gene is a member of a multi-gene family which shares strong similarity with the Drosophila dishevelled gene, dsh. The Drosophila dishevelled gene encodes a cytoplasmic phosphoprotein that regulates cell proliferation.




订购热线: 4008-898-798


Applications

Immunohistochemistry

Predicted cell location: Cytoplasm Positive control: Human prostate cancer

Recommended dilution: 25-100





The image on the left is immunohistochemistry of paraffin-embedded Human prostate cancer tissue using ml262731(DVL3 Antibody) at dilution 1/20, on the right is treated with synthetic peptide. (Original magnification: ×200)

Predicted cell location: Cytoplasm Positive control: Human brain Recommended dilution: 25-100





The image on the left is immunohistochemistry of paraffin-embedded Human brain tissue using ml262731(DVL3 Antibody) at dilution 1/20, on the right is treated with synthetic peptide. (Original magnification: ×200)

ELISA

Recommended dilution: 2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net 网址: www.mlbio.cn