

Anti-GFRA4 antibody

 Cat. No.
 ml262757

 Package
 25 μl/100 μl/200 μl

 Storage
 -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview	
Description	Anti-GFRA4 rabbit polyclonal antibody
Applications	ELISA, IHC
Immunogen	Synthetic peptide of human GFRA4
Reactivity	Human
Content	0.4 mg/ml
Host species	Rabbit
lg class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification
Target information	
Symbol	GFRA4
Full name	GDNF family receptor alpha 4

Q9GZZ7

Target Background

Synonyms Swissprot


The protein encoded by this gene is a member of the GDNF receptor family. It is a glycosylphosphatidylinositol(GPI)-linked cell surface receptor for persephin, and mediates activation of the RET tyrosine kinase receptor. This gene is a candidate gene for RET-associated diseases. Alternatively spliced transcript variants encoding different isoforms have been described for this gene.



订购热线: 4008-898-798

Applications Immunohistochemistry

Predicted cell location: Secreted or Cell membrane Positive control: Human breast cancer Recommended dilution: 10-50





The image on the left is immunohistochemistry of paraffin-embedded Human breast cancer tissue using ml262757(GFRA4 Antibody) at dilution 1/20, on the right is treated with synthetic peptide. (Original magnification: ×200)

ELISA

Recommended dilution: 1000-2000

Predicted cell location: Secreted or Cell membrane Positive control: Human esophagus cancer Recommended dilution: 10-50





The image on the left is immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using ml262757(GFRA4 Antibody) at dilution 1/20, on the right is treated with synthetic peptide. (Original magnification: ×200)

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695,2881505696、

邮箱: mlbio_cn@yeah.net 网址: www.mlbio.cn