

订购热线: 4008-898-798

# Anti-CHODL antibody

**Cat. No.** ml263452

**Package** 25 μl/100 μl/200 μl

**Storage** -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

**Product overview** 

**Description** Anti-CHODL rabbit polyclonal antibody

Applications ELISA, IHC

Immunogen Synthetic peptide of human CHODL

ReactivityHuman, MouseContent0.6 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

**Target information** 

Symbol CHODL Full name chondrolectin

**Synonyms** MT75; PRED12; C21orf68

Swissprot Q9H9P2

### **Target Background**

This gene encodes a type I membrane protein with a carbohydrate recognition domain characteristic of C-type lectins in its extracellular portion. In other proteins, this domain is involved in endocytosis of glycoproteins and exogenous sugarbearing pathogens. This protein localizes predominantly to the perinuclear region. Several transcript variants encoding a few different isoforms have been found for this gene.



订购热线: 4008-898-798


# **Applications**


## **Immunohistochemistry**

Predicted cell location: Cytoplasm and Cell membrane

Positive control: Human ovarian cancer

Recommended dilution: 10-50





The image on the left is immunohistochemistry of paraffin-embedded Human ovarian cancer tissue using ml263452(CHODL Antibody) at dilution 1/20, on the right is treated with synthetic peptide. (Original magnification: ×200)

#### **ELISA**

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio\_cn@yeah.net 网址: www.mlbio.cn