

订购热线: 4008-898-798

Anti-CITED4 antibody

Cat. No. ml263466

Package 25 μl/100 μl/200 μl

Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview

Description Anti-CITED4 rabbit polyclonal antibody

Applications ELISA, IHC

Immunogen Synthetic peptide of human CITED4

ReactivityHumanContent0.9 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

Target information

Symbol CITED4

Full name Cbp/p300 interacting transactivator with Glu/Asp rich carboxy-terminal domain 4

Synonyms

Swissprot Q96RK1

Target Background

The protein encoded by this intronless gene belongs to the CITED family of transcriptional coactivators that bind to several proteins, including CREB-binding protein (CBP) and p300, via a conserved 32 aa C-terminal motif, and regulate gene transcription. This protein also interacts with transcription factor AP2 (TFAP2), and thus may function as a co-activator for TFAP2. Hypermethylation and transcriptional downregulation of this gene has been observed in oligodendroglial tumors with deletions of chromosomal arms 1p and 19q, and associated with longer recurrence-free and overall survival of patients with oligodendroglial tumors.





订购热线: 4008-898-798

Applications

Immunohistochemistry

Predicted cell location: Cytoplasm Positive control: Human liver cancer Recommended dilution: 30-150





The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using ml263466(CITED4 Antibody) at dilution 1/65, on the right is treated with synthetic peptide. (Original magnification: ×200)

ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net 网址: www.mlbio.cn