

订购热线: 4008-898-798

Anti-CKS1B antibody

Cat. No. ml263468

Package 25 μl/100 μl/200 μl

Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview

Description Anti-CKS1B rabbit polyclonal antibody

Applications ELISA, IHC

Immunogen Synthetic peptide of human CKS1B

ReactivityHuman, MouseContent0.6 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

Target information

Symbol CKS1B

Full name CDC28 protein kinase regulatory subunit 1B

Synonyms CKS1; ckshs1; PNAS-16; PNAS-18

Swissprot P61024

Target Background

CKS1B protein binds to the catalytic subunit of the cyclin dependent kinases and is essential for their biological function. The CKS1B mRNA is found to be expressed in different patterns through the cell cycle in HeLa cells, which reflects a specialized role for the encoded protein. At least two transcript variants have been identified for this gene, and it appears that only one of them encodes a protein.




订购热线: 4008-898-798


Applications

Immunohistochemistry

Predicted cell location: Nucleus and Cytoplasm

Positive control: Human breast cancer Recommended dilution: 20-100





The image on the left is immunohistochemistry of paraffin-embedded Human breast cancer tissue using ml263468(CKS1B Antibody) at dilution 1/30, on the right is treated with synthetic peptide. (Original magnification: ×200)

ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net 网址: www.mlbio.cn