

订购热线: 4008-898-798

# Anti-ARCN1 antibody

**Cat. No.** ml263487

**Package** 25 μl/100 μl/200 μl

**Storage** -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

**Product overview** 

**Description** Anti-ARCN1 rabbit polyclonal antibody

Applications ELISA, IHC

Immunogen Synthetic peptide of human ARCN1

Reactivity Human, Mouse, Rat

Content 1.1 mg/ml Host species Rabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

**Target information** 

Symbol ARCN1
Full name archain 1
Synonyms COPD
Swissprot P48444

### **Target Background**

This gene maps in a region, which include the mixed lineage leukemia and Friend leukemia virus integration 1 genes, where multiple disease-associated chromosome translocations occur. It is an intracellular protein. Archain sequences are well conserved among eukaryotes and this protein may play a fundamental role in eukaryotic cell biology. It has similarities to heat shock proteins and clathrin-associated proteins, and may be involved in vesicle structure or trafficking.





订购热线: 4008-898-798

### **Applications**

## **Immunohistochemistry**

Predicted cell location: Cytoplasm Positive control: Human liver cancer Recommended dilution: 40-200





The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using ml263487(ARCN1 Antibody) at dilution 1/50, on the right is treated with synthetic peptide. (Original magnification: ×200)

#### **ELISA**

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio\_cn@yeah.net 网址: www.mlbio.cn