

订购热线: 4008-898-798

Anti-BRD8 antibody

Cat. No. ml263540

Package 25 μl/100 μl/200 μl

Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview

Description Anti-BRD8 rabbit polyclonal antibody

Applications ELISA, IHC

Immunogen Synthetic peptide of human BRD8

ReactivityHumanContent0.4 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

Target information

Symbol BRD8

Full name bromodomain containing 8
Synonyms SMAP; p120; SMAP2

Swissprot Q9H0E9

Target Background

The protein encoded by this gene interacts with thyroid hormone receptor in a ligand-dependent manner and enhances thyroid hormone-dependent activation from thyroid response elements. This protein contains a bromodomain and is thought to be a nuclear receptor coactivator. Multiple alternatively spliced transcript variants that encode distinct isoforms have been identified.



订购热线: 4008-898-798


Applications


Immunohistochemistry

Predicted cell location: Nucleus

Positive control: Human prostate cancer

Recommended dilution: 25-100





The image on the left is immunohistochemistry of paraffin-embedded Human prostate cancer tissue using ml263540(BRD8 Antibody) at dilution 1/20, on the right is treated with synthetic peptide. (Original magnification: ×200)

ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net 网址: www.mlbio.cn