

订购热线: 4008-898-798

# Anti-ORAOV1 antibody

**Cat. No.** ml263594

Package 25  $\mu$ l/100  $\mu$ l/200  $\mu$ l

**Storage** -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

**Product overview** 

**Description** Anti-ORAOV1 rabbit polyclonal antibody

Applications ELISA, IHC

Immunogen Synthetic peptide of human ORAOV1

ReactivityHumanContent1.1 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

**Target information** 

Symbol ORAOV1

Full name oral cancer overexpressed 1

SynonymsTAOS1SwissprotQ8WV07

### **Target Background**

Oral Cancer Overexpressed 1 (ORAOV1) is a candidate protooncogene locating on 11q13. Recent studies show that ORAOV1 acts as a primary driving force behind 11q13 gene amplification and plays a functional role in the tumorigenesis in a variety of human squamous cell carcinomas (SCCs).



订购热线: 4008-898-798

#### **Applications**

## **Immunohistochemistry**

Predicted cell location: Nucleus and Cytoplasm Positive control: Human esophagus cancer

Recommended dilution: 20-100







The image on the left is immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using ml263594(ORAOV1 Antibody) at dilution 1/20, on the right is treated with synthetic peptide. (Original magnification: ×200)

Predicted cell location: Nucleus and Cytoplasm

Positive control: Human breast cancer

Recommended dilution: 20-100





The image on the left is immunohistochemistry of paraffin-embedded Human breast cancer tissue using ml263594(ORAOV1 Antibody) at dilution 1/20, on the right is treated with synthetic peptide. (Original magnification: ×200)

#### **ELISA**

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio\_cn@yeah.net 网址: www.mlbio.cn