

Anti-DCAF5 antibody

Cat. No.	ml263679
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-DCAF5 rabbit polyclonal antibody
Applications	ELISA, IHC
Immunogen	Synthetic peptide of human DCAF5
Reactivity	Human, Mouse
Content	0.7 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	DCAF5
Full name	DDB1 and CUL4 associated factor 5
Synonyms	BCRG2; BCRP2; WDR22; D14S1461E
Swissprot	Q96JK2

Target Background

May function as a substrate receptor for CUL4-DDB1 E3 ubiquitin-protein ligase complex.

订购热线: 4008-898-798


Applications

Immunohistochemistry

Predicted cell location: Cytoplasm

Positive control: Human liver cancer

Recommended dilution: 20-100



The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using ml263679(DCAF5 Antibody) at dilution 1/20, on the right is treated with synthetic peptide. (Original magnification: $\times 200$)

ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn