

Anti-DNAJC9 antibody

Cat. No.	ml263729
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-DNAJC9 rabbit polyclonal antibody
Applications	ELISA, IHC
Immunogen	Synthetic peptide of human DNAJC9
Reactivity	Human
Content	0.54 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	DNAJC9
Full name	DnaJ heat shock protein family (Hsp40) member C9
Synonyms	JDD1; SB73; HDJC9
Swissprot	Q8WXX5

Target Background

May play a role as co-chaperone of the Hsp70 family proteins HSPA1A, HSPA1B and HSPA8.

订购热线: 4008-898-798


Applications

Immunohistochemistry

Predicted cell location: Nucleus

Positive control: Human prostate cancer

Recommended dilution: 10-50



The image on the left is immunohistochemistry of paraffin-embedded Human prostate cancer tissue using ml263729(DNAJC9 Antibody) at dilution 1/20, on the right is treated with synthetic peptide. (Original magnification: $\times 200$)

ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn