

Anti-GBGT1 antibody

Cat. No.	ml263970
Package	25 μl/100 μl/200 μl
Storage	-20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview	
Description	Anti-GBGT1 rabbit polyclonal antibody
Applications	ELISA, IHC
Immunogen	Synthetic peptide of human GBGT1
Reactivity	Human
Content	1.14 mg/ml
Host species	Rabbit
lg class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification
Target information	
Symbol	GBGT1
Full name	globoside alpha-1,3-N-acetylgalactosaminyltransferase 1 (FORS bloc

globoside alpha-1,3-N-acetylgalactosaminyltransferase 1 (FORS blood group) FS; A3GALNT; UNQ2513 Q8N5D6

Target Background

30000

Synonyms Swissprot


This gene encodes a glycosyltransferase that plays a role in the synthesis of Forssman glycolipid (FG), a member of the globoseries glycolipid family. Glycolipids such as FG form attachment sites for the binding of pathogens to cells; expression of this protein may determine host tropism to microorganisms. Alternative splicing results in multiple transcript variants.



订购热线: 4008-898-798

Applications Immunohistochemistry

Predicted cell location: Cytoplasm and Nucleus Positive control: Human thyroid cancer Recommended dilution: 30-150





The image on the left is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using ml263970(GBGT1 Antibody) at dilution 1/40, on the right is treated with synthetic peptide. (Original magnification: ×200)

ELISA

Recommended dilution: 5000-10000

Predicted cell location: Cytoplasm and Nucleus Positive control: Human breast cancer Recommended dilution: 30-150





The image on the left is immunohistochemistry of paraffin-embedded Human breast cancer tissue using ml263970(GBGT1 Antibody) at dilution 1/40, on the right is treated with synthetic peptide. (Original magnification: ×200)

- 联系电话: 4008-898-798, 021-61725725
- 联系QQ: 2881505695,2881505696、
- 邮箱: mlbio_cn@yeah.net 网址: www.mlbio.cn